
Antibodies Product Datasheet

Product Name:	PARL Monoclonal Antibody
Catalog Number:	N0504
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of PARL (aa112-167) expressed in E. Coli.
Concentration:	
Gene Name:	PARL
Protein Name:	Presenilins-associated rhomboid-like protein mitochondrial
Other Name:	PARL; PSARL; Presenilins-associated rhomboid-like protein; mitochondrial; Mitochondrial intramembrane cleaving protease PARL
Human Swiss Prot No.:	Q9H300
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PARL Monoclonal Antibody detects endogenous levels of PARL protein.
Storage and Stability:	-20°C/1 year