

## Antibodies Product Datasheet

<b>Product Name:</b>	Akt3 Monoclonal Antibody
<b>Catalog Number:</b>	N0020
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of Akt3 expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	AKT3
<b>Protein Name:</b>	RAC-gamma serine/threonine-protein kinase
<b>Other Name:</b>	AKT3; PKBG; RAC-gamma serine/threonine-protein kinase; Protein kinase Akt-3; Protein kinase B gamma; PKB gamma; RAC-PK-gamma; STK-2
<b>Human Swiss Prot No.:</b>	Q9Y243
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
<b>Specificity:</b>	Akt3 Monoclonal Antibody detects endogenous levels of Akt3 protein.
<b>Storage and Stability:</b>	-20°C/1 year