
Antibodies Product Datasheet

Product Name:	AFP Monoclonal Antibody
Catalog Number:	N0015
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human AFP expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	AFP
Protein Name:	Alpha-fetoprotein
Other Name:	AFP; HPAFP; Alpha-fetoprotein; Alpha-1-fetoprotein; Alpha-fetoglobulin
Human Swiss Prot No.:	P02771
Formulation:	Antibody are purified by protein G affinity chromatography. Liquid in PBS containing 50% glycerol and 0.03% sodium azide.
Specificity:	AFP Monoclonal Antibody detects endogenous levels of AFP protein.
Storage and Stability:	-20°C/1 year