

## Antibodies Product Datasheet

<b>Product Name:</b>	A-FABP Monoclonal Antibody
<b>Catalog Number:</b>	N0014
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of A-FABP expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	FABP4
<b>Protein Name:</b>	Fatty acid-binding protein, adipocyte
<b>Other Name:</b>	FABP4; Fatty acid-binding protein; adipocyte; Adipocyte lipid-binding protein; ALBP; Adipocyte-type fatty acid-binding protein; A-FABP; AFABP; Fatty acid-binding protein 4
<b>Human Swiss Prot No.:</b>	P15090
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
<b>Specificity:</b>	A-FABP Monoclonal Antibody detects endogenous levels of A-FABP protein.
<b>Storage and Stability:</b>	-20°C/1 year