

Antibodies Product Datasheet

Product Name:	ABCG2 Monoclonal Antibody
Catalog Number:	N0007
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ABCG2 expressed in E. Coli.
Concentration:	
Gene Name:	ABCG2
Protein Name:	ATP-binding cassette sub-family G member 2
Other Name:	ABCG2; ABCP; BCRP; BCRP1; MXR; ATP-binding cassette sub-family G member 2; Breast cancer resistance protein; CDw338; Mitoxantrone resistance-associated protein; Placenta-specific ATP-binding cassette transporter; CD antigen CD338
Human Swiss Prot No.:	Q9UNQ0
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	ABCG2 Monoclonal Antibody detects endogenous levels of ABCG2 protein.
Storage and Stability:	-20°C/1 year