

Antibodies Product Datasheet

Product Name:	P2RX7 Antibody
Catalog Number:	C51196
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal P2RX7
Concentration:	1 mg/ml
Gene Name:	P2RX7
Protein Name:	P2X purinoceptor 7 (P2X7) (ATP receptor) (P2Z receptor) (Purinergic receptor)
Other Name:	P2X purinoceptor 7 (P2X7) (ATP receptor) (P2Z receptor) (Purinergic receptor)
Human Swiss Prot No.:	Q99572
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of P2RX7 at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year