
Antibodies Product Datasheet

Product Name:	OTC Antibody
Catalog Number:	C45357
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human OTC
Concentration:	1 mg/ml
Gene Name:	OTC
Protein Name:	OTC
Other Name:	Ornithine carbamoyltransferase, mitochondrial (EC 2.1.3.3) (Ornithine transcarbamylase) (OTCase)
Human Swiss Prot No.:	P00480
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.60% sodium azide.
Specificity:	This antibody detects endogenous levels of OTC at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year