

---

## Antibodies Product Datasheet

<b>Product Name:</b>	NSE2 Antibody
<b>Catalog Number:</b>	C45083
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from part region of human protein
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	NSMCE2 C8orf36 MMS21
<b>Protein Name:</b>	E3 SUMO-protein ligase NSE2 (EC 6.3.2.-) (MMS21 homolog) (hMMS21) (Non-structural maintenance of chromosomes element 2 homolog) (Non-SMC element 2 homolog)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q96MF7
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	NSE2 Polyclonal Antibody detects endogenous levels of protein.
<b>Storage and Stability:</b>	-20°C/1 year