
Antibodies Product Datasheet

Product Name:	MET22 Antibody
Catalog Number:	C44735
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MET22
Concentration:	1 mg/ml
Gene Name:	METTTL22 C16orf68 LP8272
Protein Name:	MET22
Other Name:	
Human Swiss Prot No.:	Q9BUU2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MET22 at Human/Mouse
Storage and Stability:	-20°C/1 year