
Antibodies Product Datasheet

Product Name:	VIP Antibody
Catalog Number:	R12-2412
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VIP. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	VIP
Protein Name:	VIP peptides
Other Name:	VIP; VIP peptides
Human Swiss Prot No.:	P01282
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VIP Polyclonal Antibody detects endogenous levels of VIP protein.
Storage and Stability:	-20°C/1 year