

Antibodies Product Datasheet

Product Name:	RFC3 Antibody
Catalog Number:	R12-2332
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RFC3. AA range:178-227
Concentration:	1 mg/ml
Gene Name:	RFC3
Protein Name:	Replication factor C subunit 3
Other Name:	RFC3; Replication factor C subunit 3; Activator 1 38 kDa subunit; A1 38 kDa subunit; Activator 1 subunit 3; Replication factor C 38 kDa subunit; RF-C 38 kDa subunit; RFC38
Human Swiss Prot No.:	P40938
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RFC3 Polyclonal Antibody detects endogenous levels of RFC3 protein.
Storage and Stability:	-20°C/1 year