

Antibodies Product Datasheet

Product Name:	CBLB Antibody
Catalog Number:	C42886
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CBLB RNF56 Nbla00127
Protein Name:	E3 ubiquitin-protein ligase CBL-B (EC 6.3.2.-) (Casitas B-lineage lymphoma proto-oncogene b) (RING finger protein 56) (SH3-binding protein CBL-B) (Signal transduction protein CBL-B)
Other Name:	
Human Swiss Prot No.:	Q13191
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CBLB Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year