

Antibodies Product Datasheet

Product Name:	CALCOCO2 Antibody
Catalog Number:	C42844
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CALCOCO2
Concentration:	1 mg/ml
Gene Name:	CALCOCO2 NDP52
Protein Name:	CALCOCO2
Other Name:	Calcium-binding and coiled-coil domain-containing protein 2 (Antigen nuclear dot 52 kDa protein) (Nuclear domain 10 protein NDP52) (Nuclear domain 10 protein 52) (Nuclear dot protein 52)
Human Swiss Prot No.:	Q13137
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CALCOCO2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year