

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | AOC3 Antibody |
| Catalog Number: | C42595 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | WB, ELISA |
| Dilutions: | WB 1:500-2000 ELISA 1:5000-20000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human protein . at AA range: 510-590 |
| Concentration: | 1 mg/ml |
| Gene Name: | AOC3 VAP1 |
| Protein Name: | Membrane primary amine oxidase (EC 1.4.3.21) (Copper amine oxidase) (HPAO) (Semicarbazide-sensitive amine oxidase) (SSAO) (Vascular adhesion protein 1) (VAP-1) |
| Other Name: | |
| Human Swiss Prot No.: | Q16853 |
| Formulation: | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Specificity: | AOC3 Polyclonal Antibody detects endogenous levels of protein. |
| Storage and Stability: | -20°C/1 year |