

Antibodies Product Datasheet

Product Name:	NFIA Antibody
Catalog Number:	C42324
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NFIA
Concentration:	1 mg/ml
Gene Name:	NFIA KIAA1439
Protein Name:	Nuclear factor 1 A-type (NF1-A) (Nuclear factor 1/A) (CCAAT-box-binding transcription factor) (CTF) (Nuclear factor I/A) (NF-I/A) (NFI-A) (TGGCA-binding protein)
Other Name:	
Human Swiss Prot No.:	Q12857
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human NFIA
Storage and Stability:	-20°C/1 year