

Antibodies Product Datasheet

Product Name:	MRP4 Antibody
Catalog Number:	C42306
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MRP4
Concentration:	1 mg/ml
Gene Name:	ABCC4 MRP4
Protein Name:	Multidrug resistance-associated protein 4 (ATP-binding cassette sub-family C member 4) (MRP/cMOAT-related ABC transporter) (Multi-specific organic anion transporter B) (MOAT-B)
Other Name:	
Human Swiss Prot No.:	O15439
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MRP4
Storage and Stability:	-20°C/1 year