
Antibodies Product Datasheet

Product Name:	HA-Tag Antibody
Catalog Number:	C42226
Reactivity:	Species independent
Applications:	WB
Dilutions:	WB: 1:5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of HA-Tag
Concentration:	
Gene Name:	HA Tag
Protein Name:	HA-Tag
Other Name:	YM3202
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects C-terminal, internal, and N-terminal HA-tag fusion proteins.
Storage and Stability:	-20°C/1 year