

## Antibodies Product Datasheet

<b>Product Name:</b>	ATG5 Rabbit Antibody
<b>Catalog Number:</b>	C42085
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section)
<b>Dilutions:</b>	IHC-p 1:50-300
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of ATG5
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ATG5 APG5L ASP
<b>Protein Name:</b>	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)
<b>Other Name:</b>	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein); YM3780
<b>Human Swiss Prot No.:</b>	Q9H1Y0
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	The antibody detects endogenous ATG5 protein
<b>Storage and Stability:</b>	-20°C/1 year