

Antibodies Product Datasheet

Product Name:	DNMT1 Antibody
Catalog Number:	C42013
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DNMT1 AIM CXXC9 DNMT
Protein Name:	DNA (cytosine-5)-methyltransferase 1 (Dnmt1) (EC 2.1.1.37) (CXXC-type zinc finger protein 9) (DNA methyltransferase Hsal) (DNA MTase Hsal) (M.Hsal) (MCMT)
Other Name:	
Human Swiss Prot No.:	P26358
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DNMT1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year