
Antibodies Product Datasheet

Product Name:	ZO1 Antibody
Catalog Number:	C41994
Reactivity:	Human, Mouse
Applications:	WB, IF, ELISA
Dilutions:	WB 1:500-2000 IF 1:100-300 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	TJP1 ZO1
Protein Name:	Tight junction protein ZO-1 (Tight junction protein 1) (Zona occludens protein 1) (Zonula occludens protein 1)
Other Name:	
Human Swiss Prot No.:	Q07157
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ZO1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year