

Antibodies Product Datasheet

Product Name:	P-glycoprotein 1 Antibody
Catalog Number:	R12-2304
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human P-glycoprotein 1. AA range:534-583
Concentration:	1 mg/ml
Gene Name:	ABCB1
Protein Name:	Multidrug resistance protein 1
Other Name:	p-pg;ABCB1; MDR1; PGY1; Multidrug resistance protein 1; ATP-binding cassette sub-family B member 1; P-glycoprotein 1; CD antigen CD243
Human Swiss Prot No.:	P08183
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	P-glycoprotein 1 Polyclonal Antibody detects endogenous levels of P-glycoprotein 1 protein.
Storage and Stability:	-20°C/1 year