

Antibodies Product Datasheet

Product Name:	Akt3 Antibody
Catalog Number:	C40144
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of RAC-gamma serine/threonine-protein kinase
Concentration:	
Gene Name:	AKT3
Protein Name:	RAC-gamma serine/threonine-protein kinase
Other Name:	AKT3; PKBG; RAC-gamma serine/threonine-protein kinase; Protein kinase Akt-3; Protein kinase B gamma; PKB gamma; RAC-PK-gamma; STK-2; YM3251
Human Swiss Prot No.:	Q9Y243
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Akt3 protein.
Storage and Stability:	-20°C/1 year