

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | LC3B Antibody |
| Catalog Number: | C40132 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IF/ICC, IHC-p |
| Dilutions: | WB: 1:1000-2000 IHC: 1:200-500 IF 1:200 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Recombinant Protein of MAP LC3 β |
| Concentration: | |
| Gene Name: | MAP1LC3B |
| Protein Name: | Microtubule-associated proteins 1A/1B light chain 3B |
| Other Name: | MAP1LC3B; MAP1ALC3; Microtubule-associated proteins 1A/1B light chain 3B; Autophagy-related protein LC3 B; Autophagy-related ubiquitin-like modifier LC3 B; MAP1 light chain 3-like protein 2; MAP1A/MAP1B light chain 3 B; MAP1A/MAP1B LC3 B; Microtubule-asso |
| Human Swiss Prot No.: | Q9GZQ8 |
| Formulation: | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Specificity: | The antibody detects endogenous MAP LC3 β protein. |
| Storage and Stability: | -20°C/1 year |