
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | GST-Tag Antibody |
| Catalog Number: | C40015 |
| Reactivity: | Species independent |
| Applications: | WB |
| Dilutions: | WB: 1:5000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen. |
| Immunogen: | Recombinant Protein of GST-Tag |
| Concentration: | |
| Gene Name: | |
| Protein Name: | |
| Other Name: | YM3205 |
| Human Swiss Prot No.: | |
| Formulation: | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Specificity: | The antibody detects C-terminal, internal, and N-terminal GST fusion proteins. |
| Storage and Stability: | -20°C/1 year |