

Antibodies Product Datasheet

Product Name:	HA-tag Antibody
Catalog Number:	C40008
Reactivity:	Species independent
Applications:	WB, IP, ELISA
Dilutions:	Western Blot: 1/1000 - 1/3000. ELISA: 1/20000. IP: 1:300 Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	HA synthetic peptide conjugated to KLH.
Concentration:	1 mg/ml
Gene Name:	HA Tag
Protein Name:	HA-Tag
Other Name:	HA epitope tag; hemagglutinin
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.05% sodium azide.
Specificity:	HA-tag Polyclonal Antibody detects HA-tagged recombinant proteins or HA-tagged proteins overexpressed in cells.
Storage and Stability:	-20°C/1 year