

Antibodies Product Datasheet

Product Name:	ACO2 Antibody
Catalog Number:	C31065
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ACO2. at AA range: 421-470
Concentration:	1 mg/ml
Gene Name:	ACO2
Protein Name:	ACO2
Other Name:	Aconitate hydratase, mitochondrial (Aconitase) (EC 4.2.1.3) (Citrate hydro-lyase)
Human Swiss Prot No.:	Q99798
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACO2 Polyclonal Antibody detects endogenous levels of ACO2
Storage and Stability:	-20°C/1 year