

Antibodies Product Datasheet

Product Name:	IL-19 Antibody
Catalog Number:	C30616
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000,IHC-p 1:500-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 51-100
Concentration:	1 mg/ml
Gene Name:	IL19 ZMDA1
Protein Name:	Interleukin-19 (IL-19) (Melanoma differentiation-associated protein-like protein) (NG.1)
Other Name:	Interleukin-19 (IL-19) (Melanoma differentiation-associated protein-like protein) (NG.1)
Human Swiss Prot No.:	Q9UHD0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous IL-19
Storage and Stability:	-20°C/1 year