
Antibodies Product Datasheet

Product Name:	CD223 Antibody
Catalog Number:	C30549
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LAG3. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	LAG3
Protein Name:	Lymphocyte activation gene 3 protein
Other Name:	LAG3; FDC; Lymphocyte activation gene 3 protein; LAG-3; Protein FDC; CD223
Human Swiss Prot No.:	P18627
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD223 Polyclonal Antibody detects endogenous levels of CD223 protein.
Storage and Stability:	-20°C/1 year