

Antibodies Product Datasheet

Product Name:	Flt3 Antibody
Catalog Number:	C30429
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FLT3. AA range:451-500
Concentration:	1 mg/ml
Gene Name:	FLT3
Protein Name:	Receptor-type tyrosine-protein kinase FLT3
Other Name:	FLT3; CD135; FLK2; STK1; Receptor-type tyrosine-protein kinase FLT3; FL cytokine receptor; Fetal liver kinase-2; FLK-2; Fms-like tyrosine kinase 3; FLT-3; Stem cell tyrosine kinase 1; STK-1; CD135
Human Swiss Prot No.:	P36888
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Flt3 Polyclonal Antibody detects endogenous levels of Flt3 protein.
Storage and Stability:	-20°C/1 year