

Antibodies Product Datasheet

Product Name:	CD63 Antibody
Catalog Number:	C30372
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD63. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	CD63
Protein Name:	CD63 antigen
Other Name:	CD63; MLA1; TSPAN30; CD63 antigen; Granulophysin; Lysosomal-associated membrane protein 3; LAMP-3; Melanoma-associated antigen ME491; OMA81H; Ocular melanoma-associated antigen; Tetraspanin-30; Tspan-30; CD63
Human Swiss Prot No.:	P08962
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD63 Polyclonal Antibody detects endogenous levels of CD63 protein.
Storage and Stability:	-20°C/1 year