
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | TFPI Antibody |
| Catalog Number: | C30304 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from the C-terminal region of human TFPI. AA range:251-300 |
| Concentration: | 1 mg/ml |
| Gene Name: | TFPI |
| Protein Name: | Tissue factor pathway inhibitor |
| Other Name: | TFPI; LACI; TFPI1; Tissue factor pathway inhibitor; TFPI; Extrinsic pathway inhibitor; EPI; Lipoprotein-associated coagulation inhibitor; LACI |
| Human Swiss Prot No.: | P10646 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | TFPI Polyclonal Antibody detects endogenous levels of TFPI protein. |
| Storage and Stability: | -20°C/1 year |