

Antibodies Product Datasheet

Product Name:	SDHB Antibody
Catalog Number:	C30270
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SDHB. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	SDHB
Protein Name:	Succinate dehydrogenase [ubiquinone] iron-sulfur subunit mitochondrial
Other Name:	SDHB; SDH; SDH1; Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial; Iron-sulfur subunit of complex II; Ip
Human Swiss Prot No.:	P21912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SDHB Polyclonal Antibody detects endogenous levels of SDHB protein.
Storage and Stability:	-20°C/1 year