

Antibodies Product Datasheet

Product Name:	PI 3-Kinase p85 β Antibody
Catalog Number:	C30256
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 1-60
Concentration:	1 mg/ml
Gene Name:	PIK3R2
Protein Name:	Phosphatidylinositol 3-kinase regulatory subunit beta (PI3-kinase regulatory subunit beta) (PI3K regulatory subunit beta) (PtdIns-3-kinase regulatory subunit beta) (Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta) (PI3-kinase subunit p85-beta)
Other Name:	Phosphatidylinositol 3-kinase regulatory subunit beta (PI3-kinase regulatory subunit beta) (PI3K regulatory subunit beta) (PtdIns-3-kinase regulatory subunit beta) (Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta) (PI3-kinase subunit p85-beta)
Human Swiss Prot No.:	O00459
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous PI 3-Kinase p85 β
Storage and Stability:	-20°C/1 year