

Antibodies Product Datasheet

Product Name:	ERp57 Antibody
Catalog Number:	C30255
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PDIA3. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	PDIA3
Protein Name:	Protein disulfide-isomerase A3
Other Name:	PDIA3; ERP57; ERP60; GRP58; Protein disulfide-isomerase A3; 58 kDa glucose-regulated protein; 58 kDa microsomal protein; p58; Disulfide isomerase ER-60;Endoplasmic reticulum resident protein 57; ER protein 57; ERp57; Endoplasmic reticulum resident protein
Human Swiss Prot No.:	P30101
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERp57 Polyclonal Antibody detects endogenous levels of ERp57 protein.
Storage and Stability:	-20°C/1 year