

Antibodies Product Datasheet

Product Name:	IKKy Antibody
Catalog Number:	C30213
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000,IHC-p 1:500-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 60-120
Concentration:	1 mg/ml
Gene Name:	IKBKG FIP3 NEMO
Protein Name:	NF-kappa-B essential modulator (NEMO) (FIP-3) (I κ B kinase-associated protein 1) (IKKAP1) (Inhibitor of nuclear factor kappa-B kinase subunit gamma) (I-kappa-B kinase subunit gamma) (IKK-gamma) (IKKG) (I κ B kinase subunit gamma) (NF-kappa-B essential modifi
Other Name:	NF-kappa-B essential modulator (NEMO) (FIP-3) (I κ B kinase-associated protein 1) (IKKAP1) (Inhibitor of nuclear factor kappa-B kinase subunit gamma) (I-kappa-B kinase subunit gamma) (IKK-gamma) (IKKG) (I κ B kinase subunit gamma) (NF-kappa-B essential modifi
Human Swiss Prot No.:	Q9Y6K9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous IKKy
Storage and Stability:	-20°C/1 year