

Antibodies Product Datasheet

Product Name:	Angptl4 Antibody
Catalog Number:	C30090
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ANGPTL4. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	ANGPTL4
Protein Name:	Angiopoietin-related protein 4
Other Name:	ANGPTL4; ARP4; HFARP; PGAR; PP1158; PSEC0166; Angiopoietin-related protein 4; Angiopoietin-like protein 4; Hepatic fibrinogen/angiopoietin-related protein; HFARP
Human Swiss Prot No.:	Q9BY76
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Angptl4 Polyclonal Antibody detects endogenous levels of Angptl4 protein.
Storage and Stability:	-20°C/1 year