
Antibodies Product Datasheet

Product Name:	AR Antibody
Catalog Number:	D12046N
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human AR.
Concentration:	1 mg/ml
Gene Name:	AR
Protein Name:	Androgen receptor
Other Name:	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Human Swiss Prot No.:	P10275
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AR Polyclonal Antibody detects endogenous levels of AR protein.
Storage and Stability:	-20°C/1 year