

Antibodies Product Datasheet

Product Name:	TRAF2 Antibody
Catalog Number:	C21095
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human TRAF2.
Concentration:	1 mg/ml
Gene Name:	TRAF2
Protein Name:	TNF receptor-associated factor 2
Other Name:	TRAF2; TRAF3; TNF receptor-associated factor 2; E3 ubiquitin-protein ligase TRAF2; Tumor necrosis factor 2 receptor-associated protein 3
Human Swiss Prot No.:	Q12933
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TRAF2 Polyclonal Antibody detects endogenous levels of TRAF2 protein.
Storage and Stability:	-20°C/1 year