
Antibodies Product Datasheet

Product Name:	USP15 Antibody
Catalog Number:	C19361
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human USP15. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	USP15
Protein Name:	Ubiquitin carboxyl-terminal hydrolase 15
Other Name:	USP15; KIAA0529; Ubiquitin carboxyl-terminal hydrolase 15; Deubiquitinating enzyme 15; Ubiquitin thioesterase 15; Ubiquitin-specific-processing protease 15; Unph-2; Unph4
Human Swiss Prot No.:	Q9Y4E8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	USP15 Polyclonal Antibody detects endogenous levels of USP15 protein.
Storage and Stability:	-20°C/1 year